
# X3™ for Triathlon

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## References

1. Wang, A., Manley, M., Serkeian, P., "Wear and Structural Fatigue Simulation of Crosslinked Ultra-High Molecular Weight Polyethylene for Hip and Knee Bearing Applications," Crosslinked and Thermally Treated Ultra-High Molecular Weight Polyethylene for Joint Replacements, ASTM STP1445, Kurtz, S.M., Gsell, R., and Martell, J., Eds., ASTM International, West Conshohocken, PA 2003, pp. 151-168.
2. Yau, S.S., Wang, A., Essner, A., Manley, M., Dumbleton, J., "Sequential Irradiation and Annealing of highly Crosslinked Polyethylenes: Resist Oxidation without Sacrificing Physical/Mechanical Properties," Transactions of the 51st Annual Meeting of the Orthopaedic Research Society, Washington, DC 2005:1670.
3. Essner, A., Yau, S.S., Schmidig, G., Wang, A., Dumbleton, J., Manley, M., Serekian, P., "Reducing Hip Wear Without Compromising Mechanical Strength: A Next Generation Crosslinked and Annealed Polyethylene." Transactions of the 5th Combined Meeting of the Orthopaedic Research Societies of Canada, USA, Japan and Europe, Calgary, Canada 2004: 80.
4. Essner, A., Herrera, L., Yau, S.S., Wang, A., Dumbleton, J.H., Manley, M.T., "Sequentially Crosslinked and Annealed UHMWPE CR Knee Wear," Transactions of the 51st Annual Meeting of the Orthopaedic Research Society, Washington, DC 2005: 312.
5. Essner, A., Zeng, H., Yau, S.S., Wang, A., Dumbleton, J.H., Manley, M.T., "Stair Climbing Knee Wear Simulation: Sequentially Crosslinked and Annealed UHMWPE," Transactions of the 30th Annual Meeting of the Society for Biomaterials, Memphis 2005: 199.
6. X3™ UHMWPE maintains mechanical properties for Tensile Yield Strength and Ultimate Tensile Strength of N2Vac™ gamma sterilized UHMWPE as measured by ASTM D638. Tensile Yield Strength was 23.2 +/- 0.4 MPa and 23.5 +/- 0.3 MPa for N2Vac™ UHMWPE and X3 UHMWPE, respectively. Ultimate Tensile Strength was 54.8 +/- 2.5 MPa and 56.7 +/- 2.1 MPa for N2Vac™ UHMWPE and X3™ UHMWPE, respectively.
7. Stryker Orthopaedics Scorpio® CR tibial inserts made from X3™ UHMWPE, 72-22-0708, show a 79% reduction in volumetric wear rate versus the same insert fabricated from N2Vac™ gamma sterilized UHMWPE, 72-2-0708. The insert tested was size 7, 8mm thick. Testing was conducted under multiaxial knee simulator (multi-station MTS knee joint simulator [1]) for five million cycles using appropriate size CoCr counterfaces, a specific type of diluted calf serum lubricant and the motion and loading conditions, representing normal walking, outlined in ISO/DIS 14243-3. Volumetric wear rates were  $34.6 \pm 1.5 \text{ mm}^3/106 \text{ cycles}$  for test samples. Test inserts were exposed to a gas plasma sterilization process. In vitro knee wear simulator tests have not been shown to quantitatively predict clinical wear performance.  
[1] Essner, A., Wang, A., Stark, C., and Dumbleton, J. H., "A Simulator for the Evaluation of Total Knee Replacement Wear," 5th World Biomaterials Congress, Toronto, Canada, May 1996, pg. 580.
8. X3™ UHMWPE maintains mechanical properties after accelerated oxidative aging. No statistical difference was found for Tensile Yield Strength, Ultimate Tensile Strength and Elongation as measured per ASTM D638 before and after exposure to ASTM F2003 accelerated aging (5 Atmospheres (ATM) of oxygen at 70°C for 14 days). Tensile Yield Strength was 23.5 +/- 0.3 MPa and 23.6 +/- 0.2 MPa, Ultimate Tensile Strength was 56.7 +/- 2.1 MPa and 56.3 +/- 2.3 MPa and Elongation was 267 +/- 7% and 266 +/- 9% before and after accelerated oxidative aging, respectively.
9. X3™ UHMWPE resists the effects of oxidation. No statistical difference was found for Tensile Yield Strength, Elongation, Crystallinity and Density as measured per ASTM D638, D3417 and D1505 before and after ASTM F2003 accelerated aging (5 ATM of oxygen at 70°C for 14 days). Tensile Yield Strength was 23.5 +/- 0.3 MPa, and 23.6 +/- 0.2 MPa Ultimate Tensile Strength was 56.7 +/- 2.1 MPa and 56.3 +/- 2.3 MPa, Elongation was 267 +/- 7% and 266 +/- 9%, Crystallinity was 61.7 +/- 0.6% and 61.0 +/- 0.5% and Density was 939.2 +/- 0.1 kg/m<sup>3</sup> before and after accelerated oxidative aging, respectively.
10. Wang, A., Essner, A., Stark, C., Dumbleton, J.H., A Biaxial Line-Contact Wear Machine for the Evaluation of Implant Materials for Total Knee Joint Replacement. Wear 225-229. 1999 701-707.
11. Stryker Orthopaedics Test Report: RD-04-116.

## Reconstructive

Hips  
Knees  
Trauma & Extremities  
Foot & Ankle  
Joint Preservation  
Orthobiologics & Biosurgery

## MedSurg

Power Tools & Surgical Accessories  
Computer Assisted Surgery  
Endoscopic Surgical Solutions  
Integrated Communications  
Beds, Stretchers & EMS  
Reprocessing & Remanufacturing

## Neurotechnology & Spine

Cranio-maxillofacial  
Interventional Spine  
Neurosurgical, Spine & ENT  
Neurovascular  
Spinal Implants

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Literature Number: TK1-64  
YN/M1 0m 11/13

Copyright© 2013 Stryker  
Printed in Japan

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